



# LABORATORY SOLUTIONS



# LABORATORY EQUIPMENT

Solution Provider for all Laboratory needs



**Thomas**  
Scientific



## SUPPLIES

Flasks  
Tubes  
Pipets  
Filters

## CHEMICALS

Buffers  
Reagents  
Solvents  
Standards

## INSTRUMENTS

Balances  
Thermometers  
Meters  
Electrodes

## EQUIPMENT

Pipettors  
Centrifuges  
Hot Plates  
Mills

## SAFETY & APPAREL

Gloves  
Apparel  
First Aid  
Spill Control



[www.thomassci.com](http://www.thomassci.com)



**Thermo Fisher**  
S C I E N T I F I C

## LAB INSTRUMENTS & EQUIPMENTS

- ❖ Baths, Chillers, Circulators
- ❖ Biological Safety Cabinets
- ❖ Centrifuges
- ❖ Freezers, Ovens & Accessories
- ❖ Incubators, CO2 Incubator & Accessories
- ❖ Microbiological Incubators & Environmental Chamber
- ❖ Microplate Dispensers, Washers, ReaderpH, ISE, Conductivity & Dissolved Oxygen meters
- ❖ Pipettes, Pipettor & Accessories
- ❖ Shakers, Hotplates, Stirrers & Accessories
- ❖ Water Purification System



**KNF** **BUCHI**

**F** **fisher scientific**  
part of Thermo Fisher Scientific

# FACILITY CONSTRUCTION



## DESIGN & BUILD

- Site Survey
- Client brief
- Layout Drawings
- Perspectives
- Equipment Room Lists
- Utility Plans
- Construction Schedules
- Project Management



## CONSULTANCY

- Architecture
- Planning
- Preparing Tenders & Specifications





## ADVANCELAB® PRO STEEL CASEWORK

Our floor mounted steel casework is the lab furniture of the future. It is flexible and allows infinite modifications at site with a "sliding knee space" concept. Construction uses electrogalvanised steel sheets with heavy epoxy powder coating. Quality hinges, door handles and drawer pulls lend to durability and years of trouble free use. We stock both epoxy and phenolic resin table tops to ensure quick ship possibilities to our ever growing export division.

## NEW SERIES

### NEW JERSEY

American Sturdy Construction



### URBANO

European Simplicity



# LAB SAFETY EQUIPMENT

## FUME HOODS



"Fume hoods are a safety device in which chemical manipulations can be performed. Potentially dangerous fumes are contained; extracted via an exhaust system and expelled to the environment where it is diluted and no longer pose as a health risk."


- General Purpose
- Polypropylene
- Stainless Steel
- Walk In
- School
- Island Bench
- Customized Units

## SAFETY STORAGE CABINETS

Safety storage cabinets are essential in any area where chemicals are stored. They provide a safe storage solution; allowing personnel time to escape in emergencies.



## GAS DISTRIBUTION SYSTEM



**GCE druva™**  
CENTRAL GAS SYSTEMS





**ANALOX**  
Sensor Technology

## GAS ALARM

The Ax60+ is a multi-gas safety device for monitoring carbon dioxide and oxygen. The CO2 and O2 sensors are interchangeable and can be fully integrated as part of a multi-point, multi-gas detection and alarm system.

The O2NE+ is a simple to use and maintain ambient oxygen depletion monitor and sensor, ideal for monitoring oxygen levels where inert gases such as N2, Ar or He pose a risk of depleting O2 levels in ambient air.

## ACCESSORIES



- 1.Exhaust Arm 2.Mobile Exhaust Unit 3.PCR Cabinet 4.Tapware 5.Sanitary Ware  
6.Laminar Flow Cabinet 7.Laboratory Seating 8.Laboratory Cart 9.Ergolab Electrical Adjustable Table  
10.Balance Table



## Singapore

Advancelab (S) Pte Ltd

No.52 Senang Crescent, Singapore  
416619  
Tel: +65 6448 8255  
Fax: +65 6448 9833  
Email: info@advancelab.com.sg  
Website: www.advancelab.com.sg



## Malaysia

Advancelab Sdn Bhd

KL Office: Unit B-1-10, Block B, No. 2,  
Jalan PJU 1A/7A, Oasis Square, Ara  
Damansara, 47301 Petaling Jaya,  
Selangor Darul Ehsan, Malaysia.  
Tel: +603 7831 0188  
Fax: +603 78310588  
Email: info-my@advancelab-global.com  
Website: www.advancelab.com.sg

Johor Factory: No. 3388, Jalan Pekeliling  
Tanjung 27/2, Kawasan Perindustrian  
Indahpura, 81000 Kulaijaya, Johor,  
Malaysia.  
Tel: +607 660 8877  
Fax: +607 660 8866



## Myanmar

Advancelab Scientific &  
Engineering Co., Ltd

No(81/2), 7th Street, Than Thu Mar Road,  
(14/1)Quarter, South Okkalapa  
Township, 11091 Yangon, Myanmar.  
Tel: +95 (1) 572393 / +95 (9) 779753802  
Fax: +95 (1) 572393  
Email: info-mm@advancelab-global.com  
Website: www.advancelab.com.sg

## Vietnam

Advancelab Pte Ltd

83/10 Bach Dang Street, Ward 2,  
Tan Binh District,  
Ho Chi Minh City 700000, Vietnam  
Tel: +84 8 2244 6300  
Fax: +84 8 3845 1974  
Email: info@advancelab.com.sg  
Website: www.advancelab.com.sg



## Indonesia

PT. Advancelab Saintifik

Jl. Boulevard Raya Blok QJ 5 No. 23  
Kelapa Gading, Jakarta Utara, 14240  
Indonesia  
Tel: +6221 45854569 / +6221 45854570  
Fax: +6221 45854570  
Email: info-id@advancelab-global.com  
Website: www.advancelab.com.sg

## Thailand

Advancelab (Thailand)  
Co., Ltd.

147 / 170 - 171 Baromrajchonnane Road,  
Arunamarin Bangkoknoi, Bangkok 10700,  
Thailand.  
Tel: +66 2 884 9480  
Fax: +66 2 884 6441  
Email: primasci@advancelab-global.com  
Website: www.advancelab.com.sg

## U.A.E (Dubai)

Advancelab FZCO

No.FZJOAB1310  
Jafza One, Jebel Ali Dubai  
United Arab Emirates  
Email: info@advancelab.com.sg  
Website: www.advancelab.com.sg